Cryopreservation of Human Mesenchymal Stem Cells for Clinical Applications: **Current Methods and Challenges**

Kar Wey Yong^{1,2} Wan Kamarul Zaman Wan Safwani¹, Feng Xu,^{2,3} Wan Abu Bakar Wan Abas,¹ Jane Ru Choi,^{1,2} and Belinda Pingguan-Murphy¹

Mesenchymal stem cells (MSCs) hold many advantages over embryonic stem cells (ESCs) and other somatic cells in clinical applications. MSCs are multipotent cells with strong immunosuppressive properties. They can be harvested from various locations in the human body (e.g., bone marrow and adipose tissues). Cryopreservation represents an efficient method for the preservation and pooling of MSCs, to obtain the cell counts required for clinical applications, such as cell-based therapies and regenerative medicine. Upon cryopreservation, it is important to preserve MSCs functional properties including immunomodulatory properties and multilineage differentiation ability. Further, a biosafety evaluation of cryopreserved MSCs is essential prior to their clinical applications. However, the existing cryopreservation methods for MSCs are associated with notable limitations, leading to a need for new or improved methods to be established for a more efficient application of cryopreserved MSCs in stem cell-based therapies. We review the important parameters for cryopreservation of MSCs and the existing cryopreservation methods for MSCs. Further, we also discuss the challenges to be addressed in order to preserve MSCs effectively for clinical applications.

Introduction

S TEM CELLS ARE IDEAL CANDIDATES for many biomedical applications, particularly cell-based therapies and regenerative medicine.^{1,2} Stem cells are divided into two broad types: embryonic stem cells (ESCs), obtained from the inner cell mass of blastocysts, and adult stem cells, particularly MSCs, found in adult tissues.^{3,4} MSCs hold great potential in clinical applications by avoiding the ethical concerns related to ESCs.^{5,6} MSCs can be isolated from various locations other than bone marrow,^{7,8} including fat,^{9,10} periosteum,¹¹ amniotic fluid,¹² and umbilical cord blood.¹³

MSCs are also capable of differentiating to multiple specific types of cells such those in bone or cartilage,^{14,15} and implantation of MSCs may overcome the problems of the organ shortage crisis. MSCs also have strong immunoregulatory and immunosuppressive effects, which support the use of MSCs in allogeneic transplantation.^{16,17} Implantation of allogeneic MSCs can avoid the graft versus host diseases (GvHD), which are often observed in recipients after allogeneic organ transplantation.18,19

To achieve the large numbers of MSCs required for clinical applications, the cells either need to be expanded in culture, or MSCs from multiple donors must be pooled. In either case, a method for long-term storage of the cells is required for off-the-shelf availability.²⁰ Cryopreservation is currently the only method to preserve cells, including MSCs, for any considerable period. Cryopreservation maintains cell functional properties and allows pooling of cells to reach the cell numbers required for clinical application.²¹ Without cryopreservation, cells have to be continuously subcultured, which may accumulate genetic changes, resulting in heterogeneity or tumorigenicity.^{22,2}

At -196°C (i.e., in liquid nitrogen), cells have no metabolic demands, thus avoiding biological variation due to genetic drift even when stored for multiple years. Further, cells are protected from infection while cryopreserved.²⁴ Therefore, cryopreservation provides cells with their specific genetic characteristics and intact function at specific passages for clinical and research purposes.²³ The modern technique of cryopreservation allows the long-term storage of living cells and tissues that offer a great potential for clinical applications, including bone marrow transplantation,²⁵ bone grafts,²⁶ blood transfusion,²⁷ bone marrow transplantation,²⁵ and *in vitro* fertilization.²⁸ Although preclinical studies have shown that implantation of cryopreserved MSCs is capable of restoring

¹Department of Biomedical Engineering, Faculty of Engineering, University of Malaya, Kuala Lumpur, Malaysia. ²Bioinspired Engineering and Biomechanics Center (BEBC), and ³The Key Library of Biomedical Information Engineering of Ministry of Education, School of Life Science and Technology, Xi'an Jiaotong University, Xi'an, P.R. China.

myocardial function²⁹ and treat intestinal inflammation,³⁰ the use of cryopreserved MSCs in clinical applications is still not well established.

Various methods have been developed for cryopreservation of MSCs, such as slow freezing and vitrification. However, these methods have their limitations. First, MSCs preserved with these methods still have the risk of cryo-injury even though cryoprotective agents (CPAs) are involved in protecting the cells against cryo-injury.³¹ Second, CPAs such as DMSO may induce undesired differentiation of MSCs to neuron-like cells.³² To address these issues, optimization of CPA use or the development of new cryopreservation methods is required.

Most recently, with advances in the micro/nano technologies, it is possible to improve the efficiency of cryopreservation while overcoming the limitations of existing cryopreservation methods. Freezing cells encapsulated in nanoliter droplets on highly hydrophobic nano-rough surfaces and the development of nontoxic nanoscale bio-inspired CPAs are the promising advances in cryopreservation.³³ But for now, MSCs can only be routinely cryopreserved with the existing methods.

There are some existing reviews on cryopreservation of stem cells.^{34–37} However, this review focuses not only on the important parameters for MSCs cryopreservation, but also discusses the limitations of the existing cryopreservation methods and the challenges to be addressed in order to preserve MSCs effectively for clinical applications.

Important Parameters for Cryopreservation of MSCs

Compared to cryopreservation of specialized cells such as osteoblasts and chondrocytes, more parameters need to be evaluated for cryopreservation of MSCs. For instance, it is important to assess the immunomodulatory and multilineage differentiation ability of MSCs after cryopreservation on top of cell viability and phenotype, which are important in cryopreservation of the specialized cells. If the immunomodulatory and differentiation ability are compromised during the cryopreservation process, MSCs might not be able to suppress the transplant rejection and differentiate into other types of cells to aid in regeneration and repair of damaged organs.

Moreover, information regarding the biosafety of cryopreserved MSCs is still not well established. After cryopreservation, stem cells might run the risk of undergoing genetic instability that may affect cellular proliferation and differentiation, eventually leading to tumorigenesis.³⁸ It has been suggested that tumorigenicity of stem cells might be caused by CPAs during cryopreservation (e.g., DMSO has the potential to induce genetic instability³⁹), whereas the biosafety of other CPAs has not yet been well established.³⁸ Therefore, the assessment of the biosafety of cryopreserved MSCs in terms of genetic stability and tumorigenic potential is essential prior to clinical applications. Furthermore, the assessment of homing/engraftment potential of cryopreserved MSCs has been explored recently.

Immunomodulatory ability

MSCs hold great potential for use within immunomodulatory therapies (e.g., GvHD) and allogeneic transplantation.^{17,19} They have the ability to secrete immunosuppressive molecules

(e.g., interleukin-4 (IL-4) and IL-10) to inhibit the activation of immune cells (e.g., dendritic cells and lymphocytes), thus preventing autoimmune activities and inflammation.^{40,41} To evaluate the use of cryopreserved MSCs in these applications, their immunomodulatory ability is determined by the immunosuppressive effect on lymphocytes or peripheral blood mononuclear cells (PBMCs).^{42–46} It has been reported that cryopreservation does not affect the ability of MSCs to suppress the proliferation of lymphocytes or PBMCs.^{42,45,46} However, Francois et al.⁴³ and Moll et al.⁴⁴ suggested a potential impairment of cryopreserved MSC immunosuppressive phenotype directly after thawing, which is associated with the activation of heat-shock proteins (e.g., Hsp 27, Hsp 47, Hsp 56, and Hsp 90) that initiate a heat-shock stress on MSCs during the thawing process.⁴³ Moreover, the loss of viability observed in both studies also directly affects the immunomodulation capacity of MSCs.⁴³ The difference in viability observed in the respective studies may be due to the different methods of cultivation and freezing used.⁴⁵

Interestingly, the immunosuppressive effect of cryopreserved MSCs on lymphocytes or PBMCs could be restored to levels comparable to fresh MSCs after post-thaw culture for a short duration (fewer than 7 days).^{43,44} These results imply that cryopreservation negatively affects the immunomodulatory properties of MSCs in a reversible manner,⁴³ and the therapeutic effects of cryopreserved MSCs should be recovered in culture prior to clinical applications. Future studies should aim to improve the therapeutic efficacy and delivery mode of cryopreserved MSCs for efficient clinical applications.⁴⁴

Cell differentiation capacity

Stemness markers such as Nanog, Rex-1, Sox-2, and Oct-4 have been used to analyze the multipotency of stem cells.^{47,48} Studies have shown that the reduction of MSCs differentiation capacity was associated with the decreased expression level of those markers.^{48,49} Therefore, stemness markers should be evaluated as they might be affected by cryopreservation and thus compromise the differentiation ability of stem cells. Angelo et al.¹² and Janz Fde et al.⁵⁰ have shown that MSCs are capable of retaining their differentiation ability by maintaining the stemness markers (e.g., Nanog, Sox-2, and Oct-4) level.

The differentiation of MSCs into specific cell types such as adipocytes, osteocytes, and chondrocytes has been assessed by histochemical staining and gene expression studies. The presence of round lipid droplets stained by Oil Red O was observed after cryopreserved MSCs were induced to differentiate into adipogenic-like cells.^{10,12,51–53} These adipogenic-like cells expressed adipogenic markers such as lipoprotein lipase and peroxisome proliferator-activated receptor- γ (PPAR- γ).^{50,51}

To assess the osteogenic potential of cryopreserved MSCs, the presence of calcium deposition after osteogenic induction of cryopreserved MSCs was confirmed by Alizarin Red or Von Kossa staining.^{10,13,51,53,54} These osteogenic-like cells showed a high expression level of osteogenic markers such as osteocalcin (OSC) and alkaline phosphatase (ALP).^{51,54} The presence of proteoglycan after chondrogenic induction of cryopreserved MSCs was confirmed by Alcian Blue or Safranin O staining.^{10,12,51,55} These chondrogenic-like cells expressed chondrogenic markers such as collagen type II.⁵⁰

Taken together, although the stemness profiles of cryopreserved MSCs have not been evaluated in most of the studies, existing evidence has not shown a significant effect of cryopreservation on the MSC differentiation ability, especially in term of adipogenic, osteogenic, and chondrogenic differentiation. On the contrary, James et al.⁵⁶ reported that adiposederived MSCs preserved in 10% DMSO and 90% fetal bovine serum (FBS) displayed low adipogenic and osteogenic potential when compared to fresh cells, by downregulating the expression of adipogenic and osteogenic genes. Therefore, it is essential to determine and compare the effects of cryopreservation on gene expression of differentiation capacity and stemness among MSCs preserved in various CPAs, to provide insight into molecular changes that may occur following cell freezing²⁰.

Cell viability

To achieve an adequate therapeutic effect from immediately applied cryopreserved MSCs, it is very important to evaluate the cell viability. CPAs are essential to maintain the cell viability when the cells are stored at $-196^{\circ}C.^{57}$ Viability assays such as trypan blue exclusion, annexin V-propidium iodide (annexin V-PI), and live-dead cell staining (e.g., acetomethoxy derivate of calcein (calcein-AM)/ethidium bromide) assays, have been used to evaluate the viability of MSCs after cryopreservation. These assays consistently indicate that DMSO ($\leq 10\%$) gives a high cell viability of MSCs (>75%) following cryopreservation.^{7,8,50,58}

In general, cell viability is affected by different types and concentrations of CPAs used in the cryopreservation process. A comparison study on cell viability of amnionderived MSCs preserved in various types of CPAs showed that DMSO gives a higher cell viability (>80%) than glycerol (70%–80%) or nonpermeating CPAs such as trehalose (30%–40%) and sucrose (20%–30%),⁵⁰ indicating the relatively high efficiency of DMSO in maintaining the survival rate of MSCs throughout the freezing and thawing process.

Cell phenotype and proliferation capacity

The cell phenotype of MSCs can be determined through their morphology and surface marker (CD marker) expression. Morphologically, MSCs have spindle or fibroblast-like shapes.^{59,60} It has been reported that there were no morphological changes to MSCs after cryopreservation.^{9,53,61} MSCs are positive for mesenchymal-associated markers such as CD90, CD105, and CD73, while they lack hematopoietic associated markers such as CD34, CD19 CD45, CD14, and HLA DRDPDQ.^{62,63} After cryopreservation, MSCs retained a similar expression level of positive and negative CD markers.^{7,50,53–55,61} Overall, existing evidence shows that cryopreservation has no significant effect on phenotype of MSCs.

The proliferation capacity of MSCs can be determined through their growth kinetics. Growth kinetics can be analyzed by calculating population doubling time.^{64,65} It has been reported that population doubling time was similar between cryopreserved and noncryopreserved MSCs from various adult tissues.^{9,55,61} These results indicate that MSCs are capable of maintaining their proliferative potential after cryopreservation.

Biosafety (genetic stability and tumorigenicity)

To determine whether stem cells are safe to be used in clinical applications, potential risks of the stem cells must be assessed through the analytical assessment of the product characteristics, animal studies and nonclinical safety assessment.⁶⁶ A failure from insufficient safety assessment resulting in tumor formation or life threatening graft rejection, indicates that they are not safe to be used in clinical applications.⁶⁷ Because genetic aberrations have been strongly associated with cancers, it is important that cell preparations for clinical use are free from cancer-associated genomic alterations, and this requires defined culture conditions and the genetic characterization of the final clinical product.^{22,68}

According to Diaferia et al.,³⁸ a cryopreservation process involving DMSO has the potential to modify the cell cycle and chromosome stability of stem cells and lead to alteration of cell functions that might eventually result in tumorigenesis. As chromosomal aberrations are hallmarks of human cancer,⁶⁹ it is very important to analyze the total choromosomal content of cells through karyotyping prior to clinical applications. Karyotyping allows the detection of aneuploidy (loss and gain of chromosomes) and chromosomal anomalies like inversions, depletions, translocation and duplications, which are frequently observed in genetic diseases.⁷⁰ The cytogenetic status of cryopreserved MSCs has been determined in previous studies using the karyotyping method, which demonstrated that cryopreservation did not alter chromosomes' numbers and structures in MSCs.^{12,45,46,51}

Tumorigenesis can also be induced through epigenetic changes to cellular function.⁷¹ Both DNA methylation and histone modification are among mechanisms that cause such changes.⁷² DMSO is likely to affect these epigenetic changes on cellular functions, particularly proliferation and differentiation of stem cells.³⁹ The effect of other CPAs on epigenetic changes has not yet been explored.³⁸ These epigenetic changes cannot be detected by the karyotyping method. Therefore, the tumorigenic potential of cryopreserved MSCs should be further evaluated through tumorigenic assessment such as DNA damage, and expression levels of both tumor suppressor genes and oncogenes. These concerns have further strengthened the need for more thorough biosafety evaluation of cryopreserved MSCs. In addition, biosafety comparison studies among various types of CPAs are needed to determine the most suitable CPA use for clinical therapy.

Homing/engraftment potential

The homing/engraftment of MSCs is the process by which cells migrate, engraft, and exert local functional effects (e.g., immunosuppresion).⁷³ For example, upon sensing an injury signal, MSCs will migrate to areas of inflammation and suppress the inflammation process.⁷⁴. However, Chinnadurai et al.⁷⁵ showed that cryopreservation reduces the homing/engraftment potential of MSCs, indicated by poor binding to the extracellular matrix (e.g., fibronectin) and human endothelial cells compared to fresh MSCs. This reduction is associated with a reduced cytoskeletal F-actin content in cryopreserved MSCs. The critical process that led to the change in F-actin content was the freezing process.⁷⁶ To date, the impact of cryopreservation on homing/engraftment potential of MSCs has not been fully evaluated, and therefore should be further explored.

Cryopreservation Methods for MSCs

There are two main methods of cryopreservation: 1) slow freezing and 2) rapid freezing/vitrification. Both methods

Та	BLE 1. COMPARISON	BETWEEN SLO	ow Freezing
	AND VITRIF	ICATION METHO	od
T -	A	Slow	Vituifogtion

No	Aspect	freezing	Vitrification
1	Concentration of CPAs*	Low	High
2	Risk of freeze injury	High	Low
3	Post-thaw viability	High	High
4	Risk of toxicity of CPAs	Low	High
5	Potential contamination with pathogenic agents	Low	High
6	Manipulation skill	Easy	Good manipulation skill is needed.

*CPAs, Cryoprotective agents.

have their advantages and disadvantages (Table 1) when it comes to MSC cryopreservation. The limitations of the existing cryopreservation methods should be addressed to ensure cryopreservation efficiency.

Apart from freezing, optimizing the thawing method of frozen MSCs is also important. To date, the standard method of thawing frozen MSCs (either from slow freezing or vitrification) is to warm them rapidly (>100°C/min) in a water bath at 37°C until all ice crystals disappear.³⁵ This method generally results in high post-thaw recovery of viable cells without using high-cost equipment.

However, due to the potential contamination of the waterbath with microorganisms, it might be safer to thaw cell suspensions using a dry warming procedure.⁷⁷ To achieve this, frozen MSCs can be thawed at 10°C/min in a controlledrate freezing/thawing chamber. This method has resulted in high post-thaw viability of MSCs comparable to those thawed with the standard method.⁷⁸ The development of a low-cost, reliable, and controlled-rate dry warming device for thawing frozen MSCs is essential for efficient clinical applications.

After thawing, MSCs are washed by centrifugation prior to clinical applications to remove the CPAs, especially those which are toxic (e.g., DMSO). However, this method has resulted in cell loss that might affect the clinical outcome. Therefore, the development of methods to remove CPAs while minimizing cell loss would benefit all recipients of cryopreserved MSC transplants.³⁷

Slow freezing

The slow freezing method is preferable for MSC cryopreservation in clinics and research laboratories, due to the low risk of contamination and easier processing.³³ With a freezing rate of 1°C/min, a large number of MSCs can be easily frozen in one vial with the required CPA at low concentrations (<1.5 M).⁷⁹ Further, there is no direct contact between cells and nonsterile liquid nitrogen during the freezing process, significantly reducing the potential risk of contamination. Initially, MSCs loaded with CPAs were frozen using a nonprogrammable time freezing protocol in a -20° C freezer or freezing device called "Mr. Frosty" that produces a freezing rate of 1°C/min. To improve cryopreservation efficiency, a programmable freezing time protocol was developed using a high-cost controlled-rate freezer. However, a comparison study showed that either of these protocols can be used to freeze human MSCs without affecting their viability, phenotypes, and functional properties.⁵⁰

With an advance in technology, a programmed freezer with alternating magnetic field and electric field called the 'Cell Alive System' (CAS) has been introduced. The principle of this system is to prevent the formation of ice by vibrating the cells and water molecules during freezing with a nonthermal mechanism.⁸⁰ This can further reduce the risk of freeze injury to MSCs.⁸¹ To date, the nonprogrammable freezing time protocol is still the preferable protocol for MSC cryopreservation due to its high cryopreservation efficiency and low cost. Optimization of the use of CPA is the only remaining challenge to be addressed in order to preserve MSCs effectively using slow freezing method. Xenofree 5% DMSO was reported to be an ideal CPA to preserve MSCs for clinical applications, after this approach was accurately compared with other CPAs in terms of preserving the viability, phenotype, and functional properties of MSCs in a quantitative manner.⁸²

Limitations of slow freezing. Slow freezing involves issues related to a higher risk of freeze injury (e.g., cell death) due to the formation of intra- and extracellular ice during the freezing process.^{83,84} To address this issue, optimization of the use of CPA is very important to avoid ice crystal formation by loading the MSCs with the optimum concentration of suitable CPAs (Fig. 1). DMSO with a concentration of 10% (v/v) combined with FBS (20%–90%) (v/v) is often used to preserve MSCs.^{13,51,56,61} DMSO can penetrate cells and remove water from them, thus preventing intracellular ice formation and cell rupture.³⁴

However, DMSO is cytotoxic at temperatures beyond 4°C.⁸⁵ It is essential to reduce the toxicity by removing the DMSO from cryopreserved cells prior to clinical use, but the total removal of DMSO is complex and time consuming.³⁵ This has resulted in the need for development of an alternative cryopreservation medium or a cryopreservation medium consisting of low concentration of DMSO.⁸⁶ In general, FBS is routinely added to stabilize the cell membranes and adjust intra- and extracellular osmotic pressure, thus maintaining the survival rate of cells. However, the use of FBS in cryopreservation media should be minimized due to its potential to trigger a xenogeneic immune response or transmit pathogens to the recipient.^{85,87}

Challenges of overcoming the limitations of slow freezing. The introduction of nontoxic disaccharides such as trehalose, and hydrophilic macromolecules such as polyvinylpyrollidone as a CPA, has the potential to replace the use of DMSO.^{50,52} These CPAs cannot move across the cell membrane but are able to protect the cells from rupture by forming a viscous glassy shell around the outer surface of the cell.⁵⁷ However, they are less efficient than DMSO in terms of maintaining the survival rate of MSCs. Therefore, the alternative CPAs or CPA cocktails used to replace DMSO completely will require further investigation.

To replace the use of FBS, xeno-free cryomedium (e.g., sericin and 5% human albumin solution plus 5% DMSO)^{10,88} have been proposed. Moreover, our recent study also demonstrates that adipose-derived MSCs preserved in cryomedium containing only 5% DMSO were able to maintain their phenotype, stemness, proliferative and differentiation potential, and high cell viability in a way comparable to 10% DMSO.⁸²



These results indicate that xeno-free 5% DMSO may be an ideal CPA to preserve MSCs effectively for clinical applications. However, further investigation is needed to evaluate its effect on the biosafety of MSCs.

Vitrification

Vitrification is a process by which cell suspensions are transformed directly from the aqueous phase to a glass state, after direct exposure to liquid nitrogen.⁸⁹ Two such approaches have been used to preserve cells. One of the approaches has been termed equilibrium vitrification. This approach requires both formulation of multimolar CPA mixtures and their introduction to the cell suspensions in a stepwise fashion, to help lower chemical toxicity before the cells were flash frozen with liquid nitrogen. However, this approach has the risk of inflicting osmotic damage on the cells.³⁶ So far, only human amnion-derived MSCs have been preserved with vitrification using this approach, and the results showed that vitrification does not affect the phenotype and functional properties of MSCs.⁹⁰

An alternative approach is nonequilibrium vitrification, which employs extremely high freezing rates in conjunction with one lower concentrations of CPA mixture for preserving the cells.³⁶ In this approach, vitrification systems are divided into carrier-based systems and carrier-free systems (Fig. 2). Various types of carrier for vitrification, such as straws, quartz microcapillaries, and cryoloops have been developed to achieve a high freezing rate and to prevent the cells from cryo-injury.

Each carrier system offers a different freezing rate. Plastic straws were the initial carrier used for vitrification, but only provide a freezing rate at 2500°C/min. Quartz micro-capillaries, which provide a higher freezing rate of 250,000°C/min were developed, followed by cryoloops offering a freezing rate as high as 700,000°C/min.⁷⁹ The use of higher freezing/thawing rates permits vitrification using lower CPA concentrations, thereby reducing the potential for CPA toxicity and osmotic damage. To achieve a higher

freezing rate, microscale (nanoliter) droplet vitrification systems (carrier-free systems) that involve generation of cell encapsulating CPA droplets, followed by direct injection into liquid nitrogen, have been developed.⁷⁹ To date, MSCs have not been preserved using nonequilibrium vitrification.

Limitations of vitrification. Vitrification carrier systems generally need manual handling with great manipulation skills, leading to a low throughput process. The approach only fits well with cryopreservation of cells in small volumes such as oocytes⁹¹ but not with MSCs in large volumes. Further, vitrification requires a high concentration of CPAs (6–8 M),



FIG. 2. Carrier-based (channel-based) and carrier-free (droplet-based) methods for vitrification. In carrier-based methods, CPAs are mixed with cells in a carrier (e.g., straw and capillary). Then the cell-loaded carrier is immersed in liquid nitrogen for freezing. In carrier-free methods, cell-loaded droplets premixed with CPAs are ejected into liquid nitrogen for freezing. Direct contact of cells with liquid nitrogen is required in both methods in order to preserve cells effectively. Adapted from Xu et al.⁹⁸

FIG. 1. Mechanisms of cell

injury during cryopreservation process. An optimal concentration of suitable CPAs

should be added to cell sus-

pensions to avoid such injury. Reproduced with permission from Martin-Ibanez et al.⁸⁴ which have the potential to cause chemical toxicity and osmotic shock to cells.⁷⁹ Moreover, vitrification might also lead to potential MSCs contamination with pathogenic agents, due to the direct exposure of the cell suspensions to nonsterile liquid nitrogen. Liquid nitrogen has been reported to have the issues of potential contamination with pathogens such as Aspergillus sp. and Hepatitis B virus.^{92,93} The contamination may come from cross-contamination from infected samples in the storage tank or nonsterile liquid nitrogen itself.⁹⁴

Challenges of overcoming the limitations of vitrification. Several attempts have been carried out to improve the efficiency of vitrification, which include the stepwise addition of CPA in small volumes with a microfluidic device to minimize osmotic shock, and the development of ejector-based microscale (nanolitre) droplet generation systems for vitrification of cells at high throughput.^{95,96} The micro-droplets of cell suspensions loaded with CPA can achieve a higher freezing rate than carrier-based systems and so minimize osmotic shock. Moreover, as the carrier-free systems allow efficient vitrification of droplets in a continuous manner, throughput capability has been significantly improved to the levels required for MSC cryopreservation. However, the recycling of the cells is labor-consuming as individual cell colonies need to be selected manually.^{79,83}

To date, the risk of contamination with pathogens should be reduced for the potential clinical use of vitrified MSCs. Given the prior concerns, the development of a sterile, closed, fully automated and high throughput system that allows MSCs to be vitrified without any direct exposure to liquid nitrogen is required.³³ For future clinical applications, liquid nitrogen and ejection-based droplet generation systems must be sterilized and operated in sterile hoods. Liquid nitrogen can be sterilized by using ultraviolet radiation or sterile polyetrafluoroethylene cartridge filteration.⁹⁷

Conclusions and Future Perspectives

It is essential to have a good cryopreservation protocol in order to preserve MSCs effectively for clinical applications. The protocol must be capable of maintaining the functional properties (e.g., immunomodulatory and differentiation ability) and survival rate of MSCs. There are many issues and challenges that need to be addressed concerning existing MSCs cryopreservation protocols. First, cryopreservation protocols (e.g., types and concentrations of CPAs) should be optimized to decrease the risk of harmful effects on MSCs. For instance, the heat transfer modeling suggested by Xu et al.⁹⁸ can be applied to optimize the current protocols in order to minimize freeze injury to MSCs due to osmotic shock.

Second, MSCs should be preserved without direct exposure to liquid nitrogen, to reduce the risk of pathogenic contamination. In addition, biosafety assessments of cryopreserved MSCs are necessary to ensure the safe use of the cells prior to clinical applications. Last but not least, the development of a sterile, closed, fully automated and high throughput system for MSC cryopreservation is required for clinical applications in the near future. Overcoming these challenges would establish the standardized cryopreservation protocol of human MSCs for future clinical applications.

Acknowledgments

This study was supported by University of Malaya, High Impact Research Grant (UM.C/HIR/MOHE/ENG/44) and

Postgraduate Research Grant (PPP) (PG007-2014A) from the Ministry of Higher Education Malaysia. FX was partially supported by the Major International Joint Research Program of China (11120101002), the National Natural Science Foundation of China (11372243), the International Science and Technology Cooperation Program of China (2013DFG02930), and the National 111 Project of China (B06024).

Author Disclosure Statement

All authors declare no financial conflict of interest.

References

- Arthur A, Zannettino A, Gronthos S. The therapeutic applications of multipotential mesenchymal/stromal stem cells in skeletal tissue repair. J Cell Physiol 2009;218:237– 245.
- Baraniak PR, McDevitt TC. Stem cell paracrine actions and tissue regeneration. Regen Med 2010;5:121–143.
- Estrov Z. Stem cells and somatic cells: Reprogramming and plasticity. Clin Lymphoma Myeloma. 2009;9:S319–328.
- Richardson SM, Hoyland JA, Mobasheri R, et al. Mesenchymal stem cells in regenerative medicine: Opportunities and challenges for articular cartilage and intervertebral disc tissue engineering. J Cell Physiol 2010;222:23–32.
- Lindroos B, Suuronen R, Miettinen S. The potential of adipose stem cells in regenerative medicine. Stem Cell Rev 2011;7:269–291.
- 6. Young FE. A time for restraint. Science 2000;287:1424.
- Ginis I, Grinblat B, Shirvan MH. Evaluation of bone marrow-derived mesenchymal stem cells after cryopreservation and hypothermic storage in clinically safe medium. Tissue Eng Part C Methods 2012;18:453–463.
- Pravdyuk AI, Petrenko YA, Fuller BJ, et al. Cryopreservation of alginate encapsulated mesenchymal stromal cells. Cryobiology 2013;66:215–222.
- De Rosa A, De Francesco F, Tirino V, et al. A new method for cryopreserving adipose-derived stem cells: An attractive and suitable large-scale and long-term cell banking technology. Tissue Eng Part C Methods 2009;15:659–667.
- Minonzio G, Corazza M, Mariotta L, et al. Frozen adiposederived mesenchymal stem cells maintain high capability to grow and differentiate. Cryobiology 2014;69:211–216.
- Ferretti C, Borsari V, Falconi M, et al. Human periosteumderived stem cells for tissue engineering applications: The role of VEGF. Stem Cell Rev 2012;8:882–890.
- Angelo PC, Ferreira AC, Fonseca VD, et al. Cryopreservation does not alter karyotype, multipotency, or NANOG/SOX2 gene expression of amniotic fluid mesenchymal stem cells. Genet Mol Res 2012;11:1002–1012.
- Liu G, Ye X, Zhu Y, et al. Osteogenic differentiation of GFP-labeled human umbilical cord blood derived mesenchymal stem cells after cryopreservation. Cryobiology 2011;63:125–128.
- Callihan P, Mumaw J, Machacek DW, et al. Regulation of stem cell pluripotency and differentiation by G protein coupled receptors. Pharmacol Ther 2011;129:290–306.
- 15. Doulatov S, Daley GQ. Development. A stem cell perspective on cellular engineering. Science 2013;342:700–702.
- De Miguel MP, Fuentes-Julian S, Blazquez-Martinez A, et al. Immunosuppressive properties of mesenchymal stem cells: Advances and applications. Curr Mol Med 2012;12: 574–591.
- Nauta AJ, Fibbe WE. Immunomodulatory properties of mesenchymal stromal cells. Blood 2007;110:3499–3506.

- Kim N, Im KI, Lim JY, et al. Mesenchymal stem cells for the treatment and prevention of graft-versus-host disease: Experiments and practice. Ann Hematol 2013;92:1295– 1308.
- Aggarwal S, Pittenger MF. Human mesenchymal stem cells modulate allogeneic immune cell responses. Blood 2005; 105:1815–1822.
- Davies OG, Smith AJ, Cooper PR, et al. The effects of cryopreservation on cells isolated from adipose, bone marrow and dental pulp tissues. Cryobiology 2014;69:342– 347.
- Magalhaes R, Anil Kumar PR, Wen F, et al. The use of vitrification to preserve primary rat hepatocyte monolayer on collagen-coated poly(ethylene-terephthalate) surfaces for a hybrid liver support system. Biomaterials 2009;30:4136– 4142.
- Goldring CE, Duffy PA, Benvenisty N, et al. Assessing the safety of stem cell therapeutics. Cell Stem Cell 2011;8:618– 628.
- Stubban C, Katkov I, Loring JF, et al. Chapter 4. Cryopreservation of human embryonic stem cells. In: Schwartz JFLLWH, ed; *Human Stem Cell Manual*. Oxford Academic Press; 2007, pp. 47–55.
- Rowley SD. Hematopoietic stem cell cryopreservation: A review of current techniques. J Hematother 1992;1:233-250.
- 25. Rowley SD., Feng Z, Chen L, et al. A randomized phase III clinical trial of autologous blood stem cell transplantation comparing cryopreservation using dimethylsulfoxide vs dimethylsulfoxide with hydroxyethylstarch. Bone Marrow Transplant 2003;31:1043–1051.
- Kofron MD, Opsitnick NC, Attawia MA, et al. Cryopreservation of tissue engineered constructs for bone. J Orthop Res 2003;21:1005–1010.
- Wagner CT, Martowicz ML, Livesey SA, et al. Biochemical stabilization enhances red blood cell recovery and stability following cryopreservation. Cryobiology 2002;45:153–166.
- Anger JT, Gilbert BR, Goldstein M. Cryopreservation of sperm: Indications, methods and results. J Urol 2003;170: 1079–1084.
- Chin SP, Poey AC, Wong CY, et al. Cryopreserved mesenchymal stromal cell treatment is safe and feasible for severe dilated ischemic cardiomyopathy. Cytotherapy 2010; 12:31–37.
- Castelo-Branco MT, Soares ID, Lopes DV, et al. Intraperitoneal but not intravenous cryopreserved mesenchymal stromal cells home to the inflamed colon and ameliorate experimental colitis. PLoS One 2012;7:e33360.
- Ragoonanan V, Hubel A, Aksan A. Response of the cell membrane-cytoskeleton complex to osmotic and freeze/ thaw stresses. Cryobiology 2010;61:335–344.
- Woodbury D, Schwarz EJ, Prockop DJ, et al. Adult rat and human bone marrow stromal cells differentiate into neurons. J Neurosci Res 2000;61:364–370.
- Guven S, Demirci U. Integrating nanoscale technologies with cryogenics: A step towards improved biopreservation. Nanomedicine (Lond) 2012;7:1787–1789.
- 34. Berz D, McCormack EM, Winer ES, et al. Cryopreservation of hematopoietic stem cells. Am J Hematol 2007;82: 463–472.
- 35. Thirumala, Goebel WS., Woods EJ. Clinical grade adult stem cell banking. Organogenesis 2009;5:143–154.
- Hunt CJ. Cryopreservation of human stem cells for clinical application: A review. Transfus Med Hemother 2011;38: 107–123.

- Fleming KK, Hubel A. Cryopreservation of hematopoietic and non-hematopoietic stem cells. Transfus Apher Sci 2006;34:309–315.
- Diaferia GR, Dessi SS, Deblasio P, et al. Is stem cell chromosomes stability affected by cryopreservation conditions? Cytotechnology 2008;58:11–16.
- 39. Iwatani M, Ikegami K, Kremenska Y, et al. Dimethyl sulfoxide has an impact on epigenetic profile in mouse embryoid body. Stem Cells 2006;24:2549–2556.
- 40. Caplan AI, Correa D. The MSC: An injury drugstore. Cell Stem Cell 2011;9:11–15.
- Caplan AI, Dennis JE. Mesenchymal stem cells as trophic mediators. J Cell Biochem 2006;98:1076–1084.
- Samuelsson H, Ringden O, Lonnies H, et al. Optimizing in vitro conditions for immunomodulation and expansion of mesenchymal stromal cells. Cytotherapy 2009;11:129–136.
- 43. Francois M, Copland IB, Yuan S, et al. Cryopreserved mesenchymal stromal cells display impaired immunosuppressive properties as a result of heat-shock response and impaired interferon-gamma licensing. Cytotherapy 2012;14: 147–152.
- 44. Moll G, Alm JJ, Davies LC, et al. Do cryopreserved mesenchymal stromal cells display impaired immuno-modulatory and therapeutic properties? Stem Cells 2014; 32:2430–2442.
- 45. Luetzkendorf J, Nerger K, Hering J, et al. Cryopreservation does not alter main characteristics of Good Manufacturing Process-grade human multipotent mesenchymal stromal cells including immunomodulating potential and lack of malignant transformation. Cytotherapy 2015;17: 186–198.
- 46. de Lima Prata K, de Santis GC, Orellana MD, et al. Cryopreservation of umbilical cord mesenchymal cells in xenofree conditions. Cytotherapy 2012;14:694–700.
- 47. Choi J, Pingguan-Murphy B., Wan Abas WA, et al. Impact of low oxygen tension on stemness, proliferation and differentiation potential of human adipose-derived stem cells. Biochem Biophys Res Commun 2014;448:218–224.
- 48. Chua KH, Zaman Wan Safwani WK, Hamid AA, et al. Retropatellar fat pad-derived stem cells from older osteoarthritic patients have lesser differentiation capacity and expression of stemness genes. Cytotherapy 2014;16: 599–611.
- 49. Wan Safwani WK, Makpol S, Sathapan S, et al. The changes of stemness biomarkers expression in human adipose-derived stem cells during long-term manipulation. Biotechnol Appl Biochem 2011;58:261–270.
- Janz Fde L, Debes Ade A, Cavaglieri Rde C, et al. Evaluation of distinct freezing methods and cryoprotectants for human amniotic fluid stem cells cryopreservation. J Biomed Biotechnol 2012;2012:649353.
- 51. Miranda-Sayago JM, Fernandez-Arcas N, Benito C, et al. Evaluation of a low cost cryopreservation system on the biology of human amniotic fluid-derived mesenchymal stromal cells. Cryobiology 2012;64:160–166.
- Thirumala S, Wu X, Gimble JM, et al. Evaluation of polyvinylpyrrolidone as a cryoprotectant for adipose tissuederived adult stem cells. Tissue Eng Part C Methods 2010; 16:783–792.
- 53. Al-Saqi SH, Saliem M, Quezada HC, et al. Defined serumand xeno-free cryopreservation of mesenchymal stem cells. Cell Tissue Bank 2015;16:181–193.
- 54. Liu G, Zhou H, Li Y, et al. Evaluation of the viability and osteogenic differentiation of cryopreserved human adiposederived stem cells. Cryobiology 2008;57:18–24.

- 55. Gonda K, Shigeura T, Sato T, et al. Preserved proliferative capacity and multipotency of human adipose-derived stem cells after long-term cryopreservation. Plast Reconstr Surg 2008;121:401–410.
- James AW, Levi B, Nelson ER, et al. Deleterious effects of freezing on osteogenic differentiation of human adiposederived stromal cells in vitro and in vivo. Stem Cells Dev 2011;20:427–439.
- 57. Karlsson JO, Toner M. Long-term storage of tissues by cryopreservation: Critical issues. Biomaterials 1996;17: 243–256.
- Thirumala S, Gimble JM, Devireddy RV. Evaluation of methylcellulose and dimethyl sulfoxide as the cryoprotectants in a serum-free freezing media for cryopreservation of adipose-derived adult stem cells. Stem Cells Dev 2010; 19:513–522.
- 59. Pittenger MF, Mackay AM, Beck SC, et al. Multilineage potential of adult human mesenchymal stem cells. Science 1999;284:143–147.
- 60. Choi JR, Pingguan-Murphy B, Wan Abas WA, et al. In situ normoxia enhances survival and proliferation rate of human adipose tissue-derived stromal cells without increasing the risk of tumourigenesis. PLoS One 2015;10: e0115034.
- 61. Zhang HT, Chen H, Zhao H, et al. Neural stem cells differentiation ability of human umbilical cord mesenchymal stromal cells is not altered by cryopreservation. Neurosci Lett 2011;487:118–122.
- 62. Dominici M, Le Blanc K, Mueller I, et al. Minimal criteria for defining multipotent mesenchymal stromal cells. The International Society for Cellular Therapy position statement. Cytotherapy 2006;8:315–317.
- 63. Choi JR, Pingguan-Murphy B, Abas WABW, et al. Hypoxia promotes growth and viability of human adiposederived stem cells with increased growth factors secretion. J Asian Sci Res 2014;4:328–338.
- 64. Martinello T, Bronzini I, Maccatrozzo L, et al. Canine adipose-derived-mesenchymal stem cells do not lose stem features after a long-term cryopreservation. Res Vet Sci 2011;91:18–24.
- 65. Seo JM, Sohn MY, Suh JS, et al. Cryopreservation of amniotic fluid-derived stem cells using natural cryoprotectants and low concentrations of dimethylsulfoxide. Cryobiology 2011;62:167–173.
- Fink DW Jr. FDA regulation of stem cell-based products. Science 2009;324:1662–1663.
- 67. Sharpe ME, Morton D, Rossi A. Nonclinical safety strategies for stem cell therapies. Toxicol Appl Pharmacol 2012;262:223–231.
- 68. Laurent LC, Ulitsky I, Slavin I, et al. Dynamic changes in the copy number of pluripotency and cell proliferation genes in human ESCs and iPSCs during reprogramming and time in culture. Cell Stem Cell 2011;8:106–118.
- 69. Hanahan D, Weinberg RA. Hallmarks of cancer: The next generation. Cell 2011;144:646–674.
- Dittmar KE, Simann M, Zghoul N, et al. Quality of cell products: Authenticity, identity, genomic stability and status of differentiation. Transfus Med Hemother 2010;37: 57–64.
- Waldmann T, Schneider R. Targeting histone modifications—Epigenetics in cancer. Curr Opin Cell Biol 2013; 25:184–189.
- Nelissen EC, van Montfoort AP, Dumoulin JC, et al. Epigenetics and the placenta. Hum Reprod Update 2011;17: 397–417.

- Sackstein R. The lymphocyte homing receptors: Gatekeepers of the multistep paradigm. Curr Opin Hematol 2005;12:444–450.
- Li Y, Chen J, Chen XG, et al. Human marrow stromal cell therapy for stroke in rat: neurotrophins and functional recovery. Neurology 2002;59:514–523.
- 75. Chinnadurai R, Garcia MA, Sakurai Y, et al. Actin cytoskeletal disruption following cryopreservation alters the biodistribution of human mesenchymal stromal cells in vivo. Stem Cell Reports 2014;3:60–72.
- Xu X, Liu Y, Cui Z, et al. Effects of osmotic and cold shock on adherent human mesenchymal stem cells during cryopreservation. J Biotechnol 2012;162:224–231.
- 77. Rollig C, Babatz J, Wagner I, et al. Thawing of cryopreserved mobilized peripheral blood—Comparison between waterbath and dry warming device. Cytotherapy 2002;4: 551–555.
- 78. Thirumala S, Zvonic S, Floyd E, et al. Effect of various freezing parameters on the immediate post-thaw membrane integrity of adipose tissue derived adult stem cells. Biotechnol Prog 2005;21:1511–1524.
- 79. Zhang X, Catalano PN, Gurkan UA, et al. Emerging technologies in medical applications of minimum volume vitrification. Nanomedicine (Lond) 2011;6:1115–1129.
- Wowk B. Electric and magnetic fields in cryopreservation. Cryobiology 2012;64:301–303.
- Lee SY, Huang GW, Shiung JN, et al. Magnetic cryopreservation for dental pulp stem cells. Cells Tissues Organs 2012;196:23–33.
- Yong KW, Pingguan-Murphy B, Xu F, et al. Phenotypic and functional characterization of long-term cryopreserved human adipose-derived stem cells. Sci Rep 2015;5:9596.
- Li Y, Tan JC, Li LS. Comparison of three methods for cryopreservation of human embryonic stem cells. Fertil Steril 2010;93:999–1005.
- Martin-Ibanez R, Hovatta O, Canals JM. Cryopreservation of human pluripotent stem cells: Are we going in the right direction? In: Katkov I, ed; *Current Frontiers in Cryobiology*. InTech; 2012. pp. 139–166.
- Zambelli A, Poggi G, Da Prada G, et al. Clinical toxicity of cryopreserved circulating progenitor cells infusion. Anticancer Res 1998;18:4705–4708.
- Son JH, Heo YJ, Park MY, et al. Optimization of cryopreservation condition for hematopoietic stem cells from umbilical cord blood. Cryobiology 2010;60:287–292.
- Tuschong L, Soenen SL, Blaese RM, et al. Immune response to fetal calf serum by two adenosine deaminasedeficient patients after T cell gene therapy. Hum Gene Ther 2002;13:1605–1610.
- Verdanova M, Pytlik R, Kalbacova MH. Evaluation of sericin as a fetal bovine serum-replacing cryoprotectant during freezing of human mesenchymal stromal cells and human osteoblast-like cells. Biopreserv Biobank 2014; 12:99–105.
- Rall WF, Fahy GM. Ice-free cryopreservation of mouse embryos at -196 degrees C by vitrification. Nature 1985; 313:573-575.
- Moon JH, Lee JR, Jee BC, et al. Successful vitrification of human amnion-derived mesenchymal stem cells. Hum Reprod 2008;23:1760–1770.
- Song YS, Adler D, Xu F, et al. Vitrification and levitation of a liquid droplet on liquid nitrogen. Proc Natl Acad Sci USA 2010;107:4596–4600.
- 92. Fountain D, Ralston M, Higgins N, et al. Liquid nitrogen freezers: A potential source of microbial contamination of

hematopoietic stem cell components. Transfusion 1997; 37:585-591.

- Tedder RS, Zuckerman MA, Goldstone AH, et al. Hepatitis B transmission from contaminated cryopreservation tank. Lancet 1995;346:137–140.
- 94. AbdelHafez F, Xu J, Goldberg J, et al. Vitrification in open and closed carriers at different cell stages: assessment of embryo survival, development, DNA integrity and stability during vapor phase storage for transport. BMC Biotechnol 2011;11:29.
- 95. Samot J, Moon S, Shao L, et al. Blood banking in living droplets. PLoS One 2011;6:e17530.
- 96. Song YS, Moon S, Hulli L, et al. Microfluidics for cryopreservation. Lab Chip 2009;9:1874–1881.
- Parmegiani L, Cognigni GE, Filicori M. Ultra-violet sterilization of liquid nitrogen prior to vitrification. Hum Reprod 2009;24:2969.

 Xu F, Moon S, Zhang X, et al. Multi-scale heat and mass transfer modelling of cell and tissue cryopreservation. Philos Trans A Math Phys Eng Sci 2010;368: 561–583.

> Address correspondence to: Wan Kamarul Zaman Wan Safwani, PhD Department of Biomedical Engineering Faculty of Engineering University of Malaya 50603 Kuala Lumpur Malaysia

> > E-mail: wansafwani@um.edu.my